NSK

OPERATION MANUAL

Please read this OPERATION MANUAL carefully before use, and file for future reference.

Light Probe

isoE-LUX

CE MADE IN JAPAN

OM-Z0265E

1. User and Indications for Use

User: Qualified Professionals

Indications for Use: This medical device is intended for guiding the light supplied from a dental motor equipped with a light source to the end point through the light guide (glass fiber), thereby projecting the light into the oral cavity.

2. Precautions for handling and operation

- Please read these precautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.


Class	Degree of Risk
 CAUTION	Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.
NOTICE	General product specification information highlighted to avoid product malfunction and performance reduction.

CAUTION

- When operating the product always consider patient safety.
- Users are responsible for the operational control, maintenance and continual inspection of this product.
- Do not disassemble or alter the light probe except as recommended by NSK in this Operation Manual.
- Do not allow any impact on to the product. Do not drop the product. Impact may lead to deformation of the product or damage to the product.
- Should the light probe function abnormally during use, stop using the light probe immediately and contact your Authorized NSK Dealer.
- Do not use the following fluids to wipe, immerse or clean the product; strong / super acid water, strong acid / alkaline chemicals, chlorine-containing solutions, solvents such as benzine or thinner.
- The light probe is delivered in a non-sterile condition and must be autoclaved prior to use.
- Perform regular function and maintenance checks (Refer to "6. Periodical Maintenance Checks").
- If the light probe has not been used for a long period, check that no abnormalities are found with the product before use.
- To avoid clinical downtime it is recommended that a spare be kept on hand in case of a breakdown during treatment.
- Avoid direct eye contact with the light. Do not point the light into patient's eyes.

NOTICE

- Sometimes dark stains are found on the optic illumination point (glass fiber end) due to a characteristic of the product. This is not a defect of the product and it does not affect the light projection capability of the optic illumination point.
- The performance of this product is dependent on the light source of an optic motor. Light transmission of the product might be reduced due to maintenance or use conditions.
- Turning off the light can take a few seconds according to the connected motor.

3. Connecting and Disconnecting the Light Probe from the Motor

3-1 Connecting

- 1) Insert the light probe direct to the motor (Fig. 1). Screw the light probe until it clicks and locks into position.
- 2) Confirm that the light probe is securely connected to the motor.


Fig. 1

3-2 Disconnecting

Hold the motor and the light probe, then pull apart.

CAUTION

- Connect ONLY to E type optic motors (ISO 3964 (EN ISO 3964)).

NOTICE

- The light probe illuminates when the motor light is activated.
- The light probe is available even while the irrigation function of the connected motor is activated, but in this case water may drip when the light probe is disconnected. When using this product, disabling the irrigation function of the motor is recommended.

4. Checking the Light Probe Before Each Use

Follow the check procedure below before use. If any abnormalities are found, stop using the light probe immediately and contact your Authorized NSK Dealer.

5. Post-use Maintenance

After each patient maintain the product as follows.

CAUTION

- Immediately after a treatment, the light probe should be cleaned and sterilized. Failure to properly maintain the light probe will deteriorate the optic illumination point (glass fiber end), which will remarkably reduce the light transmittance.
- Do not use the following fluids to wipe, immerse or clean the product; strong / super acid water, strong acid / alkaline chemicals, chlorine-containing solutions, solvents such as benzene or thinner.
- Do not lubricate this product.

5-1 Preparation

- 1) Wear eye protection, a mask and gloves to avoid infection.
- 2) Disconnect the light probe from the motor.

5-2 Cleaning

■ Manual Cleaning (External Surface)

- 1) Remove debris on the light probe surface using a brush under running water (at 38 °C or less). Do NOT use a wire brush.
- 2) Wipe off moisture on the surface using an absorbent cloth.

CAUTION

- If water enters the light probe, using an air syringe, blow air into the rear of the light probe to remove the internal water.

■ Automatic Cleaning (External Surface)

Use a thermo-disinfector in accordance with ISO 15883 (EN ISO 15883).

Refer to the thermo-disinfector manual.

⚠ CAUTION

- Always place the light probe in the basket during cleaning in the thermo-disinfector. If the light probe is inserted into an adaptor while cleaning, the light probe might be accidentally released due to water pressure of the thermo-disinfector.

5-3 Cleaning of Optic Illumination Point (glass fiber end)

Wipe clean the optic illumination point (glass fiber end) using an alcohol-immersed cotton swab. Remove all debris.

⚠ CAUTION

- Do NOT use a sharp tool to clean the optic illumination point (glass fiber end). It could damage the glass and reduce the light transmission.


Fig. 2

5-4 Disinfection

■ Manual Disinfection (External Surface)

Wipe the surface with cloth soaked in disinfectant solution such as alcohol.

■ Automatic Disinfection (External Surface)

Use a thermo-disinfector in accordance with ISO 15883 (EN ISO 15883).

Refer to the thermo-disinfector manual.

5-5 Sterilization

Sterilize the product by autoclave sterilization. After each patient sterilize as noted below.

- 1) Insert into an autoclave pouch. Seal the pouch.
- 2) Autoclavable under the conditions below.

Type	Temperature	Time
Gravity-displacement Autoclave	132°C	more than 15 min.
Pre-vacuum Autoclave	134°C	more than 3 min.

- 3) The light probe should remain in the autoclave pouch until required for use.

CAUTION

- To avoid product failure, do not use a sterilizer that exceeds a cycle temperature of 138°C, including the dry cycle. In some sterilizers, the chamber temperature may exceed 138°C. Contact the sterilizer manufacturer for detailed information about cycle temperatures.
- Clean the light probe prior to autoclaving. If debris remains on the internal surface it can get stuck and cause product failure.
- Do not autoclave the product with other instruments even when it is in a pouch. This is to prevent possible discoloration and damage to the product from chemical residue on other instruments.
- Do not heat or cool the product too quickly. Rapid change in temperature could cause damage to the product.
- Keep the product in suitable atmospheric pressure, temperature, humidity, ventilation, and sunlight. The air should be free from dust, salt and sulphur.
- Do not touch the product immediately after autoclaving as it will be very hot and must remain in a sterile condition.
- Autoclave sterilization is recommended for the product. The validity of other sterilization methods is not confirmed.

NOTICE

- NSK recommends Class B sterilizers as stated in EN 13060.

6. Periodical Maintenance Checks

Every 3 months perform periodical maintenance checks, referring to the check sheet below. If any abnormalities are found, contact your Authorized NSK Dealer.

Points to check	Details
Glass Fiber	Check that the optic illumination point (glass fiber end) is NOT damaged. Also pull the tip of the light probe gently to check it is securely connected.

7. Specifications

Model	isoE-LUX
Max. Working Part Diameter	Ø3.5mm
Optic	Glass Fiber
Use Environment	Temperature: 10 - 40°C, Humidity: 30 - 75% (No Condensation)
Transportation and Store Environment	Temperature -10 - 50°C, Humidity: 10 - 85%, Pressure: 500 - 1,060hPa

8. Symbol

This product is Autoclavable up to Max.135°C.

This product can be washed via Thermo Disinfector.

Mandatory conformance mark on many products placed on the market in the European Economic Area.

Manufacturer.

Authorized representative in the European Community.

9. Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non NSK parts installed. Replacement parts are available for seven years beyond discontinuation of the model.

10. Disposing product

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose the product for you.

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NSK Europe GmbH EC REP
NSK France SAS
NSK United Kingdom Ltd
NSK America Corp
NSK Oceania Pty Ltd